

Orthopaedic



Harley Supports

Our range of supports have been used in clinical practice for many years, providing gentle and flexible support whilst stabilising the hips, pelvis and spine helping you stay safe and banishing pain. All Harley Supports are covered with a 12 month guarantee.

All SUPPORTS fitted with touch close fastening

1 Harley Power Plus Support Belt

- Ideal for heavy tasks
- Additional side pulls
- Compressive support to the low back
- Removable back protection pad
- 21cm at back narrowing to 15cm at front



| Code | Size | Waist |
|---------|---------|--------------------|
| SP24032 | Small | 56-69cm (22-27") |
| SP24046 | Medium | 71-84cm (28-33") |
| SP24047 | Large | 86-99cm (34-39") |
| SP24048 | X Large | 102-114cm (40-45") |

2 Harley Gentle Forme Support Belt

- Slim effective compressive support
- Removable back protection pad
- 20cm at back narrowing to 13cm at front



| Code | Size | Waist |
|---------|---------|--------------------|
| SP24031 | Small | 56-69cm (22-27") |
| SP24043 | Medium | 71-84cm (28-33") |
| SP24044 | Large | 86-99cm (34-39") |
| SP24045 | X Large | 102-114cm (40-45") |

3 Harley Universal Support Belt

- Stabilises hips & pelvis
- Detachable side pulls for added support
- Uniform compression with unrestricted mobility
- 16cm wide



| Code | Size | Waist |
|---------|---------|--------------------|
| SP24106 | Small | 56-69cm (22-27") |
| SP24114 | Medium | 71-84cm (28-33") |
| SP24117 | Large | 86-99cm (34-39") |
| SP24120 | X Large | 102-114cm (40-45") |

4 Harley Sacroiliac Support Belt

- Additional stabilisation of the sacrum
- Added pelvic support
- Re-establish the joints natural movement
- Additional side pulls



| Code | Size | Hip (7" below waist) |
|---------|---------|----------------------|
| SP24249 | Small | 76-86cm (30-34") |
| SP24273 | Medium | 89-102cm (35-40") |
| SP24274 | Large | 104-117cm (41-46") |
| SP24275 | X Large | 120-132cm (47-52") |

5 Harley Correcting Shoulder Support

- Corrects round shoulders
- Relieves upper back pain
- Front fastening for easy use



| Code | Size | Chest (measure under bust) |
|---------|---------|----------------------------|
| SP24108 | Small | 76-86cm (30-34") |
| SP24156 | Medium | 89-102cm (35-40") |
| SP24157 | Large | 104-117cm (41-46") |
| SP24158 | X Large | 120-132cm (47-52") |





Able2 Orthopaedic Supports

Our range of orthopaedic supports are comfortable and supportive. They are constructed from a ventilated neoprene material which provides support to the affected area. The inside of all these braces have been lined with a toweling material giving added warmth and comfort. Designed for the retail market, they are retail packed in attractive boxes with euro hooks. All these Able2 supports should be hand washed.

1 Able2 Thumb Brace

This universal thumb brace wraps around the thumb and wrist providing support to the thumb. Available as a universal size, left or right

PR79001/L Left
PR79001/R Right

2 Able2 Wrist Splint

This wrist splint provides comfortable support. It is adjusted with four hook and loop fastenings and can be used for many wrist injuries.

| Left | Right | Size | Wrist Circ |
|--------------|--------------|---------|------------|
| PR79002/L/S | PR79002/R/S | Small | 13-16½cm |
| PR79002/L/M | PR79002/R/M | Medium | 16½-20cm |
| PR79002/L/L | PR79002/R/L | Large | 20-24cm |
| PR79002/L/XL | PR79002/R/XL | X Large | 24-28cm |

3 Able2 Elbow Support

Our elbow support just slips over the elbow to provide comfortable support.

| | Size | Elbow Circ |
|------------|---------|------------|
| PR79003/S | Small | 21-24cm |
| PR79003/M | Medium | 24-26½cm |
| PR79003/L | Large | 26½-31½cm |
| PR79003/XL | X Large | 31½-36½cm |

4 Able2 Tennis Elbow Support

Our tennis elbow support is available in a universal size. It provides pressure and support to the affected area.

PR79004

5 Able2 Stabilised Knee Brace

Our knee brace is comfortable and provides support to the knee. It has an open patella and is fastened above and below the knee.

| | Size | Knee Circ |
|------------|---------|-----------|
| PR79005/S | Small | 32-36cm |
| PR79005/M | Medium | 36-39½cm |
| PR79005/L | Large | 39½-43½cm |
| PR79005/XL | X Large | 43½-47½cm |

6 Able2 Ankle Support

Our ankle support slips easily onto the ankle to provide comfort and support.

| | Size | Ankle Circ |
|------------|---------|------------|
| PR79006/S | Small | 16-21½cm |
| PR79006/M | Medium | 21½-26½cm |
| PR79006/L | Large | 26½-31½cm |
| PR79006/XL | X Large | 31½-36½cm |

1 Comfort Cool Thumb CMC Restriction Splint LATEX FREE

Unique strap surrounds the thumb for extra CMC support.

- Made of thin (1.6mm), latex-free, perforated neoprene with a terry cloth liner that helps keep skin cool and comfortable.
- Ideal for CMC joint arthritis and/or instability, tendonitis and thumb discomfort associated with repetitive motion. May be helpful for tenosynovitis of the thumb extensors/abductors.
- Light compression helps reduce thumb pain.
- Contoured strap surrounds the thumb CMC joint and first metacarpal for direct joint support. Allows full hand function to help encourage compliance.
- The thumb sleeve that supports the MP joint can be cut away as needed.
- Hand wash and air dry.
- To size, measure circumference at MP joints.

| Left | Right | | MP Circumference |
|---------|---------|-------|-------------------------|
| PR79562 | PR79563 | Small | 6½" to 7" (17 to 18 cm) |
| PR79564 | PR79565 | Med. | 7" to 8" (18 to 20 cm) |
| PR79566 | PR79567 | Large | 8" to 9½" (20 to 24 cm) |

2 Comfort Cool D-ring Wrist Splint LATEX FREE

Cool, lightweight support provides cumulative trauma relief.

Provides semi flexible support without excessive bulk or warmth. Ideal for wrist disorders such as carpal tunnel syndrome, arthritis and tendonitis. Volar aluminium support is pre-bent to position the wrist in 15° extension and can be removed for adjustments. Circumferentially designed with double layer cotton stockinette under D-ring closures for comfort and easy application. Made of ⅛" (1.6mm) thick perforated neoprene and lined with terry cloth that feels good against the skin. Hand wash and air dry. To size, measure wrist circumference.

| Left | Right | | Wrist Circumference |
|---------|---------|--------|-----------------------|
| PR52962 | PR52963 | Small | 5¾" to 6½" (15-17 cm) |
| PR52964 | PR52965 | Medium | 6¾" to 7½" (17-19 cm) |
| PR52966 | PR52967 | Large | 7¾" to 8½" (20-22 cm) |

3 Comfort Cool Wide Wrist Wrap LATEX FREE

Full coverage for the wrist and distal forearm.

The ⅛" thick perforated neoprene provides support and compression to reduce pain and discomfort. Ideal for soft tissue injuries due to repetitive motion and arthritis. Wide design provides generous overlapping around the wrist and distal forearm. Limits the final degrees of wrist motion while allowing hand function. Neoprene can be trimmed. Soft terry cloth liner feels good against the skin. Hand wash and air dry. Universal sizing.

PR79510





1 Tennis Elbow Support

LATEX FREE

This Tennis Elbow support has a high level of patient compliance, as the stretch neoprene accommodates the changing muscle bulk of the moving forearm. With a D-ring strap and adjustable pressure pad, which is easily positioned on the arm to ease the pain of Tennis Elbow. Available in three sizes.

| | | Circumference 2.5 cm (1") Below Elbow Crease |
|-----------|--------|--|
| PR/VPS/03 | Small | 23-25cm (9-10") |
| PR/VPS/05 | Medium | 28-30cm (11-12") |
| PR/VPS/07 | Large | 33-35cm (13-14") |



2 Wristlet

LATEX FREE

A wrist support made in neoprene, which retains heat and conforms well to give a light compression to soothe painful joints. Wristlet can be cut without fraying to fit swollen joints. Wristlet is ideal for Tendon Injuries, Post Fracture, Arthritis, Scaphoid injury (Wrist injury at base of thumb), Skier's Thumb, Thumb Sprain and RSI.

| Left | Right | Size | Circ 2.5cm (1") below elbow crease In cm | Circ. 2.5cm (1") below elbow crease In inches |
|-------------|-------------|---------------|--|---|
| PR/WL/14 | PR/WL/24 | Small/Medium | 14-18 | 5½"-6½" |
| PR/WL/16 | PR/WL/26 | Medium/Large | 17-22 | 6½"-8½" |
| PR/WL/LXL/L | PR/WL/LXL/R | Large/X-Large | 23-24 | 8-9½" |



2 Freedom Wrist Brace - Navy

CONTAINS LATEX

A comfortable wrist brace, with a cut away design which allows flexion at the palmar crease and opposition of the thumb to the fingers. Soft Lycra binding prevents chafing of the thumb web. This brace is ideal for broken wrists, forearm fractures, Carpal Tunnel Syndrome, Arthritis, RSI and Tendonitis.

| Left | Right | Size | Circ. In cm | Circ. In Inches |
|------------|------------|---------|-------------|-----------------|
| PR/FWBN/12 | PR/FWBN/22 | X Small | 11.5-13 | 4½"-5" |
| PR/FWBN/13 | PR/FWBN/23 | Small | 14-15 | 5½"-6" |
| PR/FWBN/15 | PR/FWBN/25 | Medium | 16.5-18 | 6½"-7" |
| PR/FWBN/17 | PR/FWBN/27 | Large | 20-21.5 | 8-8½" |
| PR/FWBN/18 | PR/FWBN/28 | X Large | 23-24 | 9-9½" |



4 Shoulder/Arm Sling

LATEX FREE

This shoulder sling can be used for broken wrist, upper arm and forearm, with its universal design the sling can be adjusted to fit all. A one-piece sling with Velcro closure and stabilising buckle.

PR/UPS Standard 5cm x 112cm (2" x 44")

1 Neck Collar

LATEX FREE

This collar is ideal to stabilise the neck after trauma, or for long term use as a night resting collar. Made in fully washable medium density closed cell foam, which will not collapse with use. Available in four sizes

| | | Depth x Length |
|---------|------------|-----------------------|
| PR/BC/1 | Small | 7.5 x 43cm (3 x 17") |
| PR/BC/2 | Medium | 9 x 47cm (3½" x 18½") |
| PR/BC/3 | Large | 9 x 51cm (3½" x 20") |
| PR/BC/4 | Large/Deep | 10 x 51cm (4 x 20") |

2 Anklet

LATEX FREE

A pull on neoprene support for light compression and ankle support. Ideal for ankle sprains, Arthritis, sports injuries and a fractured ankle bone.

| | | Shoe Size (British) | Shoe Size (Continental) |
|-----------|---------|---------------------|-------------------------|
| PR/ANK/02 | X Small | 3-4 | 36-37 |
| PR/ANK/03 | Small | 5-6 | 38-39 |
| PR/ANK/05 | Medium | 7-8 | 41-42 |
| PR/ANK/07 | Large | 9-10 | 43-44 |
| PR/ANK/08 | X Large | 11-12 | 46-47 |

3 OA Knee Sleeve

CONTAINS LATEX

An elasticated pull on knee sleeve with a soft, cotton - lined foam knee pad, to provide light support and warmth for the arthritic knee. Darts prevent the knee sleeve bunching at the back of the knee. Knee Sleeve used for unstable kneecaps and movement in the kneecaps.

| | | Knee Circumference |
|-----------|----------|--------------------|
| PR/OAS/03 | Small | 25-28cm (10"-11") |
| PR/OAS/05 | Medium | 28-31cm (11"-12") |
| PR/OAS/07 | Large | 33-36cm (13"-14") |
| PR/OAS/08 | X Large | 38-41cm (15"-16") |
| PR/OAS/09 | XX Large | 41-43cm (16"-17") |

4 Pro Tem Belt

CONTAINS LATEX

Simple but effective, the high quality elastic and two flexible posterior stays make this our most popular belt. The edges have quality Lycra binding and a back depth of 20cm (8"), cut away to 18cm (7") at the front for increased comfort when sitting. Especially helps with muscular injuries to the lower back, Lumbago, Osteoporosis, Sciatica and sports injuries.

| | | Waist Measurements |
|----------|----------|---------------------|
| PR/TB/03 | Small | 66-74cm (26"-29") |
| PR/TB/05 | Medium | 76-84cm (30"-33") |
| PR/TB/07 | Large | 86-94cm (34"-37") |
| PR/TB/08 | X Large | 96-104cm (38"-41") |
| PR/TB/09 | XX Large | 107-114cm (42"-45") |

