

Wheelchair/Scooter Accessories



All Weather Wear	18-19
------------------	-------

Bags	20
------	----

Accessories	21-24
-------------	-------

For a comprehensive range of wheelchair cushions please see pages 138-141

1 Wheely Mac

A seam sealed, fully waterproof garment offering complete rain protection. Featuring an easy to use ring pull zip. Chest 1160mm (46") Length 1770mm (70"). Available in Navy only.

PR35160

2 Wheely Mac Sleeved

A sleeved version of the Wheely Mac. Fully waterproof and seam sealed, featuring elasticated cuffs and an easy to use ring pull zip. Chest 1160mm (46") Length 1770mm (70"). Available in Navy only.

PR35161

3 Wheely Cosy

Seam sealed, and fully waterproof the Wheely Cosy features a warm polyester fur lining and an easy to use central zip. Chest 1160mm (46") Length 1060mm (42"). Available in Navy only.

PR35162

4 Wheely Poncho

A seam sealed, fully waterproof poncho, which is quick and easy to put on during sudden showers. Featuring slots to accommodate the pushing handles. Chest 1160mm (46"). Available in Navy Only.

PR35163

5 Wheelybrella

Relax in the heat of the sun or keep off the rain when out and about. The flexible arm of the Wheelybrella allows for all angles of positioning, providing complete rain and sun protection. Featuring a simple clamp to fasten the Wheelybrella to the chair. Available in Navy only.

PR35169



1



1 Scooter Control Panel Cover

Designed to protect the electronics whilst you maintain full control of the handles and brakes. Fits virtually all scooters. Colour navy.

PR35181

2



2 Scooter Cape

For complete wind and rain protection from head to toe, including your scooter. Our capes are suitable for 3 or 4 wheel scooters. Velcro slots accommodate your mirrors whilst retaining maximum waterproofing. The clear front panel enables you to clearly view the controls whilst ensuring that you can safely use your lights. Features an easy to use single zip fastening, enabling you to climb in and out of the scooter with ease.

Colour navy. Width 64cm

PR35182/S **Small - length 112cm**

PR35182/M **Medium - length 124cm**

PR35182/L **Large - length 132cm**

3



3 Scooter Storage Cover

A light weight, seam sealed, fully waterproof scooter storage cover which protects your scooter from rain and dust. With elasticated corners to secure. Available in Navy only. Length 860mm (34"), Width 430mm (17").

PR35166

4



4 Scooter Mac Lined

For head to toe protection. Specifically designed for those who have difficulty putting on conventional coats and jackets. This tartan lined mac offers warmth in the winter and keeps you cool in the summer. With a shorter back the mac can be worn over the back of the chair or tucked inside.

Colour navy. Universal chest 116cm, length 121cm.

PR35179

1 Mobility Bag

A strong, stylish and multipurpose bag, which secures around the seat of the scooter or powerchair using long, wide detachable braided straps.

Alternatively attach it around the headrest of the scooter using the frame of the rollator or the pushing handles of the wheelchair with the short, narrow straps.

- Large lined main compartment with netted side pockets
- Detachable carry strap and padded carry handles
- Ringpull zips
- Free mini bag Clips inside or outside of the mobility bag, can be secured to the side of the wheelchair or scooter or underneath the seat of your wheelchair.

Colour black. Height 42cm, width 36, depth 14cm

PR35175



2 Mini Mobility Bag

Strong, stylish and ideal for mini scooters, walkers, rollators, children's wheelchairs or as a pannier bag. It has a front pocket which can be accessed easily from the side. Colour black. Height 32cm width 25cm, depth 10cm

PR35176



3 Under Seat Wheelchair Pouch

This ingenious device provides a convenient storage compartment under the chair and is accessible from the front by the wheelchair user. Made of durable dark blue nylon. Fastens securely and easily with buckles. Size 300H x 280W x 400mmD (12H x 11W x 16"D).

PR35077



4 Scooter Bag

Strong, durable bag designed to accommodate scooters with or without headrests. Features a spacious storage compartment and a box lid for easy access.

Size: Height 400mm (16"), Width 310mm (12"), Depth 50mm (2")

PR35167



5 Crutch Walking Stick Bag

Strong and spacious bag for a wheelchair or scooter. With a zipped box lid and two side holders for your crutches or walking sticks. Colour black.

PR35173/L

Large - Height 46cm, width 32cm, depth 11cm

PR35173/S

Small - Height 33cm, width 37cm





1 Crutch Holder for Scooter

Crutches or walking sticks sit comfortably in two metal holders and secure at the top with strong Velcro fastenings. Attaches to the rear of your scooter with a universal fitting suitable for most scooters.

PR35171



2 Scooter Shelter

The perfect solution for those who have nowhere to store their scooter. Simply drive it on to the metal base plate and pull the heavy duty, fully waterproof, metal framed structure over your scooter. Secures with a chain and padlock at the base. Features two easy access flaps for battery charging.

Length 218cm (85"), Width 98cm (38½"), Height 127cm (50")

PR35170



3 Disabled Badge Holder

Designed to store and display your European style blue badge and time card. The wallet opens to display both cards, and the time card can be adjusted without the need to remove it. This wallet is ideal to protect your badges and is useful to transport them should you change vehicle. The wallet is made from leather and opens to fully display the badges behind clear pockets. The wallet is protected by polished gilt corners.

Size when closed: 173 mm. (6.75") x 146 mm. (5.75"). Size when open: 348 mm. (13.8") x 146 mm. (5.75"). Available in two colours.

PR35207/B Black

PR35207/P Purple

1 Scooter Shopper

A metal framed shopping trolley which is fully lockable and can be flat packed for easy storage. The shopping trolley has a removable bag which is strong, durable and features a cool liner. It is also supplied with a free scooter attachment which enables the shopping trolley to be mounted to virtually all scooters. Colour black. Height 59cm, width 37cm depth 25cm.

PR35172

2 Rear Metal Scooter Basket

A collapsible metal basket ideal for storing excess shopping out of the way. Attaches to the rear of your scooter with a universal fitting suitable for most scooters. Also available is a hard wearing and durable basket bag designed to specifically fit our metal rear scooter basket. Features a detachable shoulder strap. Colour black. Height 28cm, width 38cm, depth 23cm.

PR35174 Basket only

PR35174/BAG Bag only

3 Scooter Basket Liner and Cover

Keep your shopping secure and dry. Removable and fits most baskets. Colour navy.

PR35178

4 Oxygen CarryON Wheelchair Bag

Conveniently stores oxygen tanks out of the way. The lightweight Oxygen CarryON Wheelchair Bag is ideal for persons with emphysema or other respiratory conditions. Large mesh bag measures 410 x 580mm (16 x 23") and holds one or two "E" cylinder tanks. Also features an 300mm (11") deep mask/nasal cannula pouch. Machine washable.

PR35203

5 Oxygen Bag

Heavy duty bag designed to carry an oxygen canister on a wheelchair or scooter. Features a front pocket to hold an oxygen mask. Colour black. Height 59cm, circumference 38cm

PR35177





1



2



3



4



1 Wheelchair Seat Covers

100% Pure New Wool

These Real wool Wheelchair Seat Covers have the unique ability to absorb perspiration, thus preventing the skin from becoming 'clammy'. This allows the skin to breathe more effectively; very important in pressure area care. Real wool is exothermic. Unlike polyester fleeces, wool can absorb up to 33% of its own weight in moisture without feeling damp (compared with polyester's 5%) Wool will not feel cold even when damp. Real wool has a remarkable natural resilience: with the real Lamb's wool high-density 1" deep pile this distributes pressure more evenly giving immediate comfort and relief. Real wool has a gentle touch that prevents friction burns from abrasive materials. The Wheelchair Seat Covers are pre-treated to enable the real Lambs wool to be hot washable up to 161°F (71°C) and autoclaved. The Real Wool Bed Fleeces can be washed at these temperatures up to 50 times and still continue to give softness and comfort to patients.

PR20137 Wheelchair Seat Cover 46 x 41cm (18" x 16")

**PR20138 Wheelchair Seat and Back Cover
91 x 41cm (36" x 16")**

2 Extra Secure Wheelchair Belt

Additional Security and Padding.

This belt criss-crosses in the back of the wheelchair and loops around the tilt levers, providing extra security. Consists of soft comfortable foam attached to a wide, heavy duty nylon strap. Measures 120 x 560mm (4½ x 22"). Slide the buckles on the strap to adjust length. Straps measure 1020mm (40") long. Fits comfortably around hips. Machine wash and dry. Light blue with dark blue trim.

PR35057

3 ADL Lap Tray

A durable and economical lap tray solution.

The sturdy ADL Lap Tray is ideal for use in a long-term care facility. It can take vigorous use, is waterproof and easy to clean. This 6.4mm (¼") thick tray is melamine-impregnated and both top and bottom are sealed. Depressed edge prevents small objects from sliding off the tray. Features rounded corners for safety. Secures with hook and loop straps riveted to the board. Width 597mm (23½") Length 559mm (22") Cut out Length 228mm (9") Width 355mm (14")

PR35018

4 Superhead

Folding Head Rest for Conventional Wheelchairs

Superhead is the highly adaptive headrest system suitable for all standard wheelchairs.

- Fully adjustable head support
- Provides optimum comfort and is easily adjustable
- Unique folding mechanism allows system to fold with the wheelchair.
- Adjusts to fit wheelchairs from 360mm-500mm
- Durable yet lightweight, aluminium construction
- High Quality, Fire retardant fabrics for comfort and safety
- Available in two sizes

PR35300 Medium for seat width 360-440mm (14-17")

PR35301 Large for seat width 420-500mm (16.5-19.5")

1 Clip-On Glass Holder

Our most economical holder.

The Clip-On Glass holder can be attached to wheelchairs or tubular poles with a diameter of 19 to 25mm ($\frac{3}{4}$ to 1"). Steel wire holder measures 73 x 76mm ($2\frac{1}{2}$ x 3").

PR35043

2 Drink Thing

Unique holder secures a variety of container sizes.

The Drink Thing has a two-ring system that securely holds glasses, cans and bottles ranging from 8 oz. to 44 oz. (.23 to 1.3 liters).

Mounts easily onto wheelchair with plastic coated clips. Folds down when not in use. US patent.

PR35041

3 Drink-Aide

Drink-Aide consists of an insulated water bottle, a vibration resistant flexible drinking tube and a universal clamp for wheelchair mounting. Drink-Aide allows users to drink fluids with little or no upper body movement. It is one of only a few completely hands free drinking bottles of this kind available. Drink-Aide allows for more freedom and travel.

- Frees caregiver to attend to other needs
- Users no longer have to wait for a carer to quench their thirst
- May reduce the frequency of hospitalisation due to urinary tract infections caused by inadequate hydration

With Drink-Aide an individual has control over fluid intake. This is especially important to those individuals who are aware that good hydration is important for decreasing the potential for urinary tract infections, to improve skin integrity and voice quality.

PR35076

PR35076/S 1 Straw

4 Globaleather Wheelchair Gloves

These leather gloves were specially developed by a wheelchair user. They incorporate a full thumb with padding that extends to the areas necessary for comfortable wheelchair use. They are securely fitted with a hook and loop tab and have a stretch cuff for comfort and flexibility. These gloves have been upgraded for added durability.

Made from soft, durable black leather. To determine the size required, measure the hand width across the MCP knuckles.

Code	Size	Hand Width	Age Guide
PR35051	X small	70mm ($2\frac{3}{4}$ "	7-11
PR35052	Small	76mm (3")	10-15
PR35053	Medium	83mm ($3\frac{1}{4}$ "	14-17
PR35054	Large	89mm ($3\frac{1}{2}$ "	Adult
PR35055	X large	95mm ($3\frac{3}{4}$ "	Large Adult
PR35056	XX large	101mm (4")	V Large Adult

5 Globaleather Gel-Palm Wheelchair Gloves

Developed by wheelchair users to make life more comfortable for wheelchair users. Featuring Amari palms with gel padding and stretch Lycra with an absorbent towelling backing. This state of the art wheelchair glove has open fingers and a full thumb to protect from callouses and blisters. Finished in silver grey and available in 6 sizes.

Size	Size	Hand Width	Age Guide
PR35150	X small	70mm ($2\frac{3}{4}$ "	7-11
PR35151	Small	76mm (3")	10-15
PR35152	Medium	83mm ($3\frac{1}{4}$ "	14-17
PR35153	Large	89mm ($3\frac{1}{2}$ "	Adult
PR35154	X large	95mm ($3\frac{3}{4}$ "	Large Adult
PR35155	XX large	101mm (4")	V Large Adult

