

Dressing & Comfort



Dressing Aids	118-121
---------------	---------

Footwear Aids	122-123
---------------	---------

Comfort Aids	124-126
--------------	---------

Slippers	127
----------	-----

Hip Protector	128
---------------	-----

1 Bra Angel

Designed by an Occupational Therapist, to assist women who find it difficult to bring both ends of a brassiere together when dressing.

- Simple to use
- One size fits all
- Light, durable and latex free
- Can be used with most brassiere types
- Telescopic action adjusts easily for size and shortens it for storage and traveling

PR55032



2 Body Care Long Handled Combs And Brushes

Brush or comb hair more easily with a longer handle. Ideal for persons with limited arm reach or hand movement. These brushes and combs have extended handles that provide easy access to hard-to-reach places at the back of the head. Firm bristles and contour design provide better contact when brushing or combing.

Brush	Length.
PR45041 Short	300mm (11½")
PR45042 Long	360mm (14")


Comb	Length
PR45039 Short	310mm (12")
PR45040 Long	370mm (14½")



3 Body Care Nail File Holder

Rubber handle holds nail file and fits comfortably in the hand. Ideal for people who have limited grasp. The Body Care Nail File Holder has a soft rubber handle that holds the nail file safely and securely. Measures 120mm (4¾") long.

PR45030



4 Long Handled Comb & Brush

Ergonomically designed personal care products featuring anti-slip handles. Allows you to take care of your own grooming needs. These products are designed to provide maximum reach with minimum effort. The carefully designed handle fits snugly in the hand to stay in place when in use. Easy to clean as they can be fully immersed in water. Non-absorbing, so quick drying. Latex free.

PR45014 Long Handled Comb - Length 380mm (15")

PR45015 Long Handled Brush - Length 350mm (14")





1 Shoehorn And Dressing Aid

This handy device combines a dressing stick and a shoehorn. Ideal for people whose reach, dexterity or ability to bend is limited due to arthritis or back/leg injuries. Designed with a push-pull S-hook at one end, and a long-handled shoehorn at the opposite end. Hook end can be used to help put on pants, sweaters, socks and other articles of clothing, or to reach clothes hangers on high closet rods. Device is available in two lengths, 610 or 760mm (24 or 30"). Made of durable plastic with a smooth surface and rounded edges. Handle measures 32mm (1 1/4") in diameter and has a soft, built-up foam that covers its entire length.

	Length	Weight
PR55013	610mm (24")	147g (5.2 oz.)
PR55014	760mm (30")	184g (6.5 oz.)



2 Dressing Stick

Essential for anyone who has difficulty bending and a limited range of movement, or the use of only one arm or hand. A plastic coated combination hook/pusher assists in putting on shirts and jackets, pulling up slacks or skirts, or removing socks. The small hook on the opposite end aids in pulling zippers and shoelace loops. Three versions are available.

PR55030/C Collapsible: Length 700mm (27 1/2")

Each piece measures 370mm (14 1/2") when stick is collapsed

PR55030/18 Fixed length 460mm (18")

PR55030/27 Fixed length 690mm (27")



3 Dressing Sleeves


Provides assistance to a person who has difficulty dressing. Made with slide sheet fabric to ensure that clothes easily slide on. Once clothes are on it can easily be removed. Pack includes a leg and arm pouch. Washable.

PR41024/S Small

PR41024/M Medium

PR41024/L Large

PR41024/XL X Large



4 Clear Zipper Pull

Add a discreet handle to zippers.

The Clear Zipper Pull clips to zipper tabs to make grasping easier. Clear nylon loop is inconspicuous. Place thumb or finger in loop and pull for easy opening and closing. Can be removed and placed in a pocket or purse. Loop measures 170mm (6 1/2") long. No pinch is required. Sold in package of 3.

PR55018



5 Zip Grip

The Zip Grip is ideal for persons who have difficulty grasping and pulling small zipper tabs. Just insert the peg of the Zip Grip through the hole in a standard zipper, fold in half and press to snap shut. Ring diameter measures 25mm (1"). Sold in package of 6.

PR55021



1 Sock-Assist

Unique features make this the most extraordinary sock aid on the market.

The easy-to-use Sock-Assist lays flat on the floor, accommodates feet and calves of all sizes, has a nonstick surface and pulls the sock up the calf instead of releasing it at the ankle. Its innovative design allows easy donning of socks or knee-high nylons using one or two hands. Also works with elastic compression knee high socks or stockings.

It's the unusual-looking features that make this sock aid so effective:

- Two side wings keep the sock on the Sock-Assist until it is released on the calf.
- Centre tongue helps ease the sock onto the device.
- Flexible polypropylene "gives" to accommodate socks.
- Flat bottom surface prevents the Sock-Assist from rolling.
- Wide foot bed accommodates average to wide feet.
- Contoured side walls help hold the sock on the device.
- Rear flares prevent the cord from rubbing against the skin.
- Textured surface prevent skin from sticking to the sock aid.
- Cord handle can be cut or looped to adjust length.

PR55006 With Loop Handle

PR55024 With Two Cord Handles



2 Tights Aid

A one-piece flexible aid for tights, with webbing handles attached. The plastic blade is shaped to hold the tights until they can be pulled up within reach. Very lightweight, and easy to use. Handle can be adjusted if necessary.

PR55320



3 Luxury Stocking Aid

This luxury stocking device is particularly useful for putting on socks or stockings. Strong wire, tong style easy to grip, telescopic handles, and cleverly shaped bendable gutter so socks or stockings can slip over straight on to the foot. This easy and efficient, handy device has a friction aided heel pad. Handle length 381mm (15")

PR55318

4 Sock/Stocking Aid

The unique shape of the flexible Sock Aid keeps socks on the aid until the sock is pulled all the way onto the foot. The pliable trough is easy to bend. The soft plastic aid rolls into the sock and creates an opening for the foot. The loop handle measures 910mm (36") long. Requires the use of two hands.

PR55319



5 Easy-Pull Sock Aid

Foot slides easily on the nylon lining while the sock remains secure. The flexible plastic core of the Easy-Pull Sock Aid allows people with the use of just one hand to shape the trough easily. Large loop handles at the end of the strong, 810cm (32") straps are great for persons with weak grasps. Nylon lining allows foot to slide easily along the trough. Terry cloth outer covering keeps the sock or stocking securely on the sock aid. Slits in the plastic core allows for contour bending around the heel. Compact, lightweight and hand washable.

PR55005





1 Good Grips® Button Hook

Get a secure, cushioned grip on hard-to-grasp buttons. The soft, cushioned grip of the Good Grips® Button Hook makes the sometimes difficult task of buttoning clothes easy. The built-up handle is made of rubber like material and features flexible ribbing that adapts to any grip. This handle is so unique, it carries a U.S. patent. The length of the handle measures 110mm (4¼"), and the total length measures 160mm (6¼"). Diameter of handle is 38mm (1½").

PR55020



2 Compression Stocking Aid


This Compression Stocking Aid is designed to take the struggle out of putting on Compression Stockings. It is made of steel with a smooth plastic coating. Convenient handles greatly assist in pulling the stockings up. Simply stretch the stocking over the column, insert toes and pull up with a continuous motion. Stocking Aids reduce bending and increase independence whilst dressing. The inside column width is 100mm (4").

PR55321

3 Heel Guide Compression Stocking Aid

Now the easiest way to put on compression hosiery. Unique indented channel on the inside of the centre cone acts as a heel guide and keeps the foot in perfect position while effortlessly pulling on the compression hose. Extra wide cone is ideal for people with edema or larger feet. Sturdy handles, 432mm (17") high, reduce the need for bending or stretching. Handles also feature 3 grips to accommodate different heights and sizes. Easy to assemble and can be taken apart for travelling.

PR55322



1 Plastic Shoehorn With Hook

Hooked handle offers a comfortable hold for fingers and may be used to pull up socks or aid dressing. Measures 410mm (16¼").

PR55028



2 Good Grips® Shoehorn

This shoehorn offers an extended reach and a secure, cushioned grip. The built-up handle is made of a latex free, rubber-like material and features flexible ribbing that adapts to any grip. This handle is so unique, it carries a U.S. Patent. People who have limited mobility and difficulty bending over will appreciate the extra-long design. Available in 3 different lengths. The textured surface prevents the foot from sticking to the shoehorn, helping it slide easily into the shoe. Handle diameter is 38mm (1½").

PR55010 Length 460mm (18")

PR55011 Length 610mm (24")

PR55012 Length 760mm (30")



3 Shoe Remover

Remove shoes effortlessly with the Shoe Remover. This product is ideal for persons with a limited range of motion or those who have difficulty bending over. Made of sturdy plastic. Not suitable for clients who have had total hip replacements or restricted hip abduction. Colour may vary.

PR55029



4 Oedema Socks

This cotton rich sock has an extra wide, extra deep cut, perfect for people who suffer from Oedema in the feet and legs. These oversized non-compression socks are 100% elastic free and have an extra wide expansion for a comfortable fit.

- Non-binding top
- Non-constricting
- Extra wide expansion
- Anti-bacterial finish
- Hand knitted seamless toe.

Available in navy and in two sizes.

PR35180/S Size 4-8

PR35180/L Size 9-13





1 Coilers

Easy-to-use Coilers never need tying. Coilers may be tightened or loosened, but never need tying. They have no pressure points and allow natural movement. Fits up to eight pairs of eyelets. Sold in pairs.

- PR55001 Black
- PR55002 Brown
- PR55003 White
- PR55004 Royal Blue



2 Elastic Shoelaces

Turn a lace-up shoe into a slip-on shoe.

Ideal for people with limited back mobility, decreased hip flexion or limited hand dexterity. Lace and tie shoes in the normal manner just once. Shoes can be slipped on or off easily without tying and untying. Elastic shoelaces measure 610mm (24") long. Two pairs per package.

- PR55000 Black
- PR55022 Brown
- PR55023 White



3 Deluxe Elastic Shoelaces

Durable laces for active wear.

Designed for athletic shoes with five, six or seven eyelet sets. Thicker and easier to handle than other elastic shoelaces, these laces are $\frac{1}{4}$ " (6.4mm) wide. Available in two lengths. White only. Two pairs per pack.

- PR55026 686mm (27") length
- PR55027 940mm (37") length




1 Toe Crests

Gel Toe Crests provide Gel comfort for pain caused by hammer, claw, or mallet toes. Fits comfortably between your ball-of-foot and toes, and relieves pressure from the toe area. Sold in pairs

PR55039/SM Small/medium

PR55039/ML Medium/large



2 Toe Protector Tubes

These Gel Toe Protector Tubes give all-around toe protection. Cushions and protects toes against pressure and friction. One 10cm length that can be cut/trimmed with scissors to fit toes.

PR55038/SM Small/medium

PR55038/ML Medium/large



3 Toe Caps

Gel Toe Caps are so soft and comfortable, just perfect for protecting and cushioning sore toes. Sold in pairs

PR55040/SM Small/medium

PR55040/ML Medium/large



4 Toe Spreaders

The toe spreaders are extremely comfortable due to the polymer's soft conforming properties. Use them to wedge between an overlapping toe or bunion deformity. Sold in pairs.

PR55041/SM Small/medium

PR55041/ML Medium/large





Bunion Protectors

These Bunion Protectors are made from unique polymer gel offering superior pain relief. They conform to the toe area giving all day comfort and provide soft gel support designed to protect enlarged joints and bunions. Universal size. Sold in a pair

PR55036



Toe Separators

These polymer gel toe separator pads are used to relieve pressure and keep toes apart. If you have painful corns between toes or simply wish to separate the toes then this shaped pad provides comfortable protection. Sold in pairs.

PR55037/SM Small/medium

PR55037/ML Medium/large



Pro Active Insoles

The Gel Pro-active insoles are designed to give superior cushioning and support to hard working feet and joints. They feature a thermoplastic rubber gel layer with anatomically shaped pads giving high energy return and shock absorption. Their unique arch massaging pads help relieve pain while the antimicrobial top cloth can help control odour and reduce friction known to cause blisters. Sold as a pair.

PR55042/M Mens

PR55042/L Ladies



Pedi-cure for Tough Nails

A little help for those tougher toe nail and feet tasks. Smooths and buffs calluses on feet and hands. Shapes and contours toenails and fingernails. 3 snap-on buffers (fine, medium and course).

Ergonomically designed contoured handle. 2-Way variable speed control. Uses 3 x AA batteries (not included).

PR61424



1 Bed Socks

These protect your feet from the cold but do not apply pressure or restrict circulation. They can be worn in bed or just for lounging. Supplied in pairs. Wash by hand. Do not tumble dry.

PR55033/M Medium - Size 4-6

PR55033/L Large - 7-11



2 Heel Protectors

These protect sore and delicate heels and ankles. There are no raised seams around the heel area. Supplied in pairs. Universal. Wash by hand. Do not tumble dry.

PR55034



3 Elbow Protectors

These offer protection for your elbows, perfect for those prone to scuffing and sores. Supplied in pairs. Universal. Wash by hand. Do not tumble dry.

PR55035



4 Fleece Slippers

These offer complete support and protection for your feet, perfect for those prone to scuffing and sores. Supplied in pairs. Wash by hand. Do not tumble dry.

PR55419/4-6 Small - Size 4-6

PR55419/7-9 Medium - Size 7-9

PR55419/10-13 Large - Size 10-13





Dunlop Slippers

These comfortable and stylish slippers have a wide opening for swollen feet, the hook and loop fastening ensures these slippers are easy to put on and take off. These machine washable slippers have non slip hardwearing soles. The slippers are available in a variety of sizes and in male and female versions.



Bluebell



Betsy



Albert

Code	Slipper	Male/Female	Size	Colour
PR55150/BL/3	Bluebell	Female	3	Blue
PR55150/BL/4	Bluebell	Female	4	Blue
PR55150/BL/5	Bluebell	Female	5	Blue
PR55150/BL/6	Bluebell	Female	6	Blue
PR55150/BL/7	Bluebell	Female	7	Blue
PR55150/BL/8	Bluebell	Female	8	Blue
PR55150/BU/3	Bluebell	Female	3	Burgundy
PR55150/BU/4	Bluebell	Female	4	Burgundy
PR55150/BU/5	Bluebell	Female	5	Burgundy
PR55150/BU/6	Bluebell	Female	6	Burgundy
PR55150/BU/7	Bluebell	Female	7	Burgundy
PR55150/BU/8	Bluebell	Female	8	Burgundy
PR55151/BL/3	Betsy	Female	3	Blue
PR55151/BL/4	Betsy	Female	4	Blue
PR55151/BL/5	Betsy	Female	5	Blue
PR55151/BL/6	Betsy	Female	6	Blue
PR55151/BL/7	Betsy	Female	7	Blue
PR55151/BL/8	Betsy	Female	8	Blue
PR55151/BU/3	Betsy	Female	3	Burgundy
PR55151/BU/4	Betsy	Female	4	Burgundy
PR55151/BU/5	Betsy	Female	5	Burgundy
PR55151/BU/6	Betsy	Female	6	Burgundy
PR55151/BU/7	Betsy	Female	7	Burgundy
PR55151/BU/8	Betsy	Female	8	Burgundy
PR55152/BL/6	Albert	Male	6	Blue
PR55152/BL/7	Albert	Male	7	Blue
PR55152/BL/8	Albert	Male	8	Blue
PR55152/BL/9	Albert	Male	9	Blue
PR55152/BL/10	Albert	Male	10	Blue
PR55152/BL/11	Albert	Male	11	Blue
PR55152/BL/12	Albert	Male	12	Blue
PR55152/BU/6	Albert	Male	6	Burgundy
PR55152/BU/7	Albert	Male	7	Burgundy
PR55152/BU/8	Albert	Male	8	Burgundy
PR55152/BU/9	Albert	Male	9	Burgundy
PR55152/BU/10	Albert	Male	10	Burgundy
PR55152/BU/11	Albert	Male	11	Burgundy
PR55152/BU/12	Albert	Male	12	Burgundy

1 The Hipshield

Reducing the risk of hip fractures in the elderly, or for post operative care

The Hipshield can help to prevent a fall becoming an injury and is designed to absorb and deflect impact away from the hipbone. The Hipshield itself is discreet and comfortable to wear and inside the protector is a lightweight honeycomb structure that compresses on impact. The Hipshield just slips into a pair of protector pants and is unobtrusive beneath skirts or trousers.

Pants are available in two styles, male and female, in six sizes based on waist measurement. Made of 95% cotton, each pair has internal pockets that you simply 'slip' the Hip protector into before putting on the garment. Suitable for day or nightwear the Hipshield is available in 'single' or 'triple pack'.

The single pack contains one knickers/pants plus two Hip protectors.

The triple pack contains three knickers/pants plus two Hip protectors.

The perfect solution, one pair to wear, one in the wash and one spare! For sizing see chart below:



Size	Female Waist size	Single pair	Triple pack
Small	660-710mm (26-28")	PR20300	PR20300-3
Medium	710-760mm (28-30")	PR20301	PR20301-3
Large	760-800mm (30-32")	PR20302	PR20302-3
X-Large	810-860mm (32-34")	PR20303	PR20303-3
XX-Large	860-910mm (34-36")	PR20304	PR20304-3
XXX-Large	910-960mm (36-38")	PR20320	PR20320-3

Size	Male Waist size	Single pair	Triple pack
Small	760-810mm (30-32")	PR20305	PR20305-3
Medium	830-890mm (33-35")	PR20306	PR20306-3
Large	910-960mm (36-38")	PR20307	PR20307-3
X-Large	990-1040mm (39-41")	PR20308	PR20308-3
XX-Large	1060-1120mm (42-44")	PR20309	PR20309-3
XXX-Large	1110-1160mm (44-46")	PR20321	PR20321-3

Washing

The Hipshield protectors should be hand washed with warm soapy water and be left to dry. This normally takes around 10 minutes.

The Undergarments may be washed by hand or in a washing machine at 60-90°.

